

001056 Sample vial (20mL) (10pcs)

After receiving the product, please check the contents.



Contents

Sample vial (20 mL)

10 pcs

This sample vial can be used for SVC-2 and SVC-3.
Sample volume from 5 to 10 mL.

Please wash the vial thoroughly before experiment.
Any occurrence of damage after use is out of our
guarantee.

Option

002012 GCE Glassy carbon electrode OD:6mm ID:3mm

012167 RE-1B Reference electrode (Ag/AgCl)

012177 Sample holder dia 9 mm

012670 Teflon cap for SVC-2

012671 Teflon cap for SVC-3

001209 Cell holder for 20 mL vial

BAS Inc.

<https://www.als-japan.com> email: sales@als-japan.com