


011224 VC-4 Voltammetry cell

After receiving the product, please check the contents.



Contents

(A) Purge tube (ETFE)	30 cm
(B) Sample vial (5mL)	7 pcs
(C) Platinum counter electrode 5.7cm	1 pc
(D) Cell holder for 5mL vial	1 pc
(E) Teflon cap for VC-4	1 pc

You can use working and reference electrodes up to 6 mm in diameter for this kit.

Please wash the vial thoroughly before using.

If the product is damaged after use, we cannot guarantee it.

Optional products

- 002012 GCE Glassy carbon electrode OD:6mm ID:3mm
- 012167 RE-1B Reference electrode (Ag/AgCl)

BAS Inc.

<https://www.als-japan.com>

email: sales@als-japan.com