
014050 RRDE-3A/CS-3A gas purge tube repair kit


Thank you for your purchase from BAS Inc. This product is a tube repair kit for RRDE-3A Rotating Ring Disk Electrode Apparatus or CS-3A Cell Stand. Please check the detailed instructions for tube replacement on our website below before proceeding. https://www.als-japan.com/support-instrument-manual.html

Contents




	(A) Silicone tube (B) Silicone tube + Teflon tube	1 рс
/	(with black cover)	1 pc
1	(C) Silicone tube + Teflon tube	1 pc
	WDIeces turn off the neuron and unrities	

%Please turn off the power and unplug the power cable before performing replacement.



CS-3A piping tube replacement



- (A): Tube from gas inlet connector to needle valve
- (B): Purge tubes from the needle valve through the pinch valve (with black cover)
- (B): Purge tubes from the needle valve through the pinch valve

Relative Product

012779 CS-3A Cell Stand Ver.1.1

013725 RRDE-3A Rotating Ring Disk Electrode Apparatus Ver.3.0 Previous product

RRDE-3A Rotating Ring Disk Electrode Apparatus Ver.2.0

BAS Inc. https://www.als-japan.com email: sales@als-japan.com